




REMOVABLE CHEST STRAP KIT PN 427039

Installation Procedure

CONTENTS

- Webbing Straps (2)
- Female Quick Release Buckle (1)
- Male Quick Release Buckle (1)



Installation Procedure

Compatible with: Black Diamond • Libra • Black Ice

1. Locate both right and left inward facing chest strap retainers on the BC shoulder straps (*Fig. 1, Step 1*).

2. With the BC facing you, hold one webbing strap (sewn fold toward you) and insert the strap from the inside to the outside of the BC (*Fig. 1, Step 2*).

3. Insert the same strap end through its own chest strap retainer from front to back. Pull on the loose end (*Fig. 1, Step 3*).

4. Repeat steps 1 through 3 for the other side.

5. Attach the male quick release buckle to the BC's left chest strap (installed in steps 1-3), by holding the male buckle with the ribbed center bar forward and the prongs pointing left (*Fig. 2*).




Figure 1




Figure 2

6. Thread the chest strap from back to front through the inside buckle slot (*Figs. 2 & 3*).




Figure 3

7. Loop the strap around the ribbed center bar and thread from front to back through the outside buckle slot (*Fig. 4*).




Figure 4

8. Pull on the loose end to tighten. Just prior to diving, the straps can be snugged for final adjustment.

9. Repeat steps 5 through 8 to attach the female quick release buckle to the BC's right chest strap.



If you have any questions regarding the information found in this manual, please contact your regional Aqua Lung Dealer or Distributor. Distributor information is available on the Aqua Lung website at: www.aqualung.com

